

TOROBlue® RT Kit for qPCR-FD-

【Catalogue Number】 RTQ-108

【Packing Information】 120 reactions for a total 20 μ L reaction volume
(or 240 reactions for a total 10 μ L reaction volume)

【Description】

TOROBlue® RT Kit for qPCR-FD- is a room-temperature stable RT kit to perform for qPCR. The engineered reverse transcriptase in this kit allows the synthesis of high-quality first-strand cDNA. 2 \times RT buffer for qPCR is optimized for highly efficient synthesis of short-chain cDNAs suitable for qPCR. The protocol is simple, and the reaction can be completed in 10-25 min.

【Feature】

- Room-Temperature stable:** the performance is not easily decrease during during storing and shipping.
- Easy-to-use:** add 50 μ L 2 \times RT buffer into a single tube to prepare the RT premix easily.
- Maximum flexibility:** cutting each tube allows for 5 or 10 RT reactions.
- High performance:** an engineered reverse transcriptase allows the synthesis of high-quality cDNA.

【Components】

TOROBlue® RT Kit for qPCR-FD- includes the following reagents, which can be used for 120 reactions for a total 20 μ L reaction volume or 240 reactions for a total 10 μ L reaction volume .

NO.	Components	Size
1	Freeze-dried RT PreMix-002-	8-tube strips \times 3
2	2 \times RT buffer for qPCR	1.25mL/tube \times 1
3	Water-DEPC treated	1.25mL/tube \times 1

Notes:

-Freeze-dried RT PreMix-002 contains lyophilized Reverse transcriptase, RNase Inhibitor and , Random hexamer and/or oligo(dT)20 primers; 2 \times RT buffer for qPCR contains reaction buffer optimized for highly efficient synthesis of short-chain cDNAs suitable for qPCR.

【Protocol】

1. Preparation of the 2 \times RT premix

-Open the reagent packaging and check if the freeze-dried powder in the tube is deliquescent.

Notes: It is not recommended to use moisture absorbing test tube.

-Cut open the 8-strips tube according to the required quantity.

Notes: The remaining test tubes should still be sealed in a dry environment.

-Add 50 μ L 2 \times RT buffer into each tube to prepare the 2 \times RT premix.

-Mix gently with a pipette and briefly centrifuged.

2. Preparation of RNA template.

-Use the purified RNA templates after quality assessment and genomic-DNA contamination assessment according to the MIQE guidelines.

-Instead of using 2 \times RT premix, 2 \times RT buffer as no-reverse transcription control was used to assessed the extent of genomic-DNA contamination in the following steps.

Notes: If gDNA contamination affects the Cq values, it is essential to eliminate it by DNase treatment.

-RNA diluent should be prepared in a thin-walled PCR tube as follows.

Component	Volume
Water-DEPC treated	10-X μ L
Total RNA(up to 3 μ g)	X μ L
Total	10 μ L

3. Preparation of the reaction mix.

-Prepare the following reaction mix in the above reaction tube with RNA diluent.

Component	Volume
RNA diluent (Step 2)	10 μ L
2 \times RT premix (Step 1)	10 μ L
Total	20 μ L

-Mix gently with a pipette and briefly centrifuged.

4. RT reaction Condition

37°C, 15min

50°C, 5min

98°C, 2 min

-Store the reaction tube one ice.

-Transfer 100 μ L RNase free water into the above reaction tube.

-Gently vortexed and briefly centrifuged.

-The cDNA products can be used directly or after dilution for real-time PCR.

【Storage】

Store the test tube at 2-8°C in a dry environment for 24 months.

【References】

Bustin SA, Benes V, Garson JA, etc,al. The MIQE guidelines: minimum information for publication of quantitative real-time PCR experiments.ClinChem.2009,Apr;55(4):611-22.

【Application】

- [1] Xu K , Meng Z , Xian X M ,et al.LncRNA PVT1 induces chondrocyte apoptosis through upregulation of TNF- α in synoviocytes by sponging miR-211-3p[J].Molecular and Cellular Probes, 2020, 52:101560.
- [2] Su X , Zhu Y , Hou Y ,et al.PVT1 induces NSCLC cell migration and invasion by regulating IL-6 via sponging miR-760.[J].Molecular and cellular probes, 2020, 54:101652.
- [3] Xia CP , Pan T , Zhang N,et al.Sp1 promotes dental pulp stem cell osteoblastic differentiation through regulating noggin[J]. Molecular and cellular probes,2020, 50:101504.
- [4] MY Xie,LH Yang,JY Cheng,et al.Gracillin relieves pulmonary fibrosis by suppressing the STAT3 axis[J]. Journal of Ethnopharmacology,2023,316:116704,
- [5] Wang Y, Fan J, Tong Y ,et al. Bioinformatics analysis of ferroptosis-related gene AKR1C3 as a potential biomarker of asthma and its identification in BEAS-2B cells. Comput Biol Med. 2023 ,3(158):106740.
- [6] Wen J , Jiang F , Liu M ,et al.Identification and expression analysis of Cathepsin B-like protease 2 genes in tomato at abiotic stresses especially at High temperature[J].Scientia Horticulturae, 2021, 277(10):109799.

【Contact】



TOROIVD TECHNOLOGY COMPANY LIMITED.

Head Office: Building #C1, 880 JiangYang Nan Road, Baoshan,Shanghai,China. 200439.

Factory: Building #20, 888 Zhujiang Road, Rudong ,Jiangshu,China. 226400

Tel: +86-21-68030217

Mail: market@toroivd.com

http: //www.toroivd.com